



Exclusive Distributor

Benetec Advanced Medical Systems

Kelderbeemd 1A, 2470 Retie

Email: info@benetecmed.com

Tel: +32 (0)14 37 79 20

SINGLE USE FLEXIBLE CYSTOSCOPE



**SINGLE USE STILL
NEEDS PERFORMANCE**

BRIGHT & CLEAR

HIGH IMAGE QUALITY

INNOVEX Single-Use Flexible Cystoscope provides every patient a clean scope with high quality view - every time. Eliminate the cost and risk associated with reprocessing and cross-contamination.



Dual deflection
210°up / 180°down

Metal snake bone structure
Excellent maneuverability
Working length 380mm

Photo & video button

Working channel luer port

Lever lock

Latest CMOS technology
Superior image
Field of view 110°±10%
Depth of field: 5mm-100mm

110°±10%

4 LED illumination

5.5mm/2.4mm

CLEAR & SHARP IMAGE

DMR-UB1-013-01 A/0